**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**F +49 921 912-780**  
**ortho@medi.de**  
**www.medi.de**

medi Australia Pty Ltd  
 Unit 4/13-21 Hallmark Street  
 Pendle Hill NSW 2145  
 Australia  
 T +61 2 9890 8696  
 F +61 2 9890 8439  
 sales@medi.au  
 www.medi.australia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: +32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 104-1375 Lionel-Boulet,  
 Varennes, Québec,  
 QC Canada J3X 1P7  
 T +1 450-583-3317 / +1 800 361 3153  
 F +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

**MAXIS a.s.,**  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

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
medi Trading (Shanghai) Co. Ltd.  
 Room 613, HnR Group Building  
 No.1178, Tian Yao Qiao Road  
 Xuhui District, Shanghai, 200232  
 People's Republic of China  
 T +86-21 50582318  
 infocn@medi.cn  
 www.medi.cn

Medi Turk Ortopedi Medikal Ic ve Dis Tic. AS  
 Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5  
 06510 Çankaya Ankara  
 Turkey  
 T: +90 312 435 20 26  
 F: +90 312 434 17 67  
 info@medi-turk.com  
 www.medi-turk.com

UK Responsible Person  
 medi UK Ltd.  
 Plough Lane  
 Hereford HR4 0EL  
 Great Britain  
 T +44-1432 37 35 00  
 F +44-1432 37 35 10  
 enquiries@mediuk.co.uk  
 www.mediuk.co.uk

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#### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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Polyester, Polyamide, Polyurethan, Federstahl

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# medi Smartfit OA

## Intended purpose

medi Smartfit OA is a knee brace for relief  and guidance  of the knee joint.

## Indications

All indications that require load relief of a compartment via the 3-point pressure relief principle, additional physiological support and, where required, limitation of the range of motion of the knee joint, such as:

- Unilateral osteoarthritis
- Cartilage surgery or chondroplasty
- Meniscus intervention

## Contra-indications

Combined lateral and medial osteoarthritis

## Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paresthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material. The immobilisation which is achievable by the device presents a risk for development of thrombosis in the lower veins of the legs, especially in the context of surgical interventions or trauma.

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the

intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Setting the restrictions for flexion and extension (can only be carried out by an orthopaedic technician)

Extension restriction wedges

(stretch restriction): -5°, 0°, 10°, 20°, 30°, 45°

Flexion restriction wedges (bend restriction) 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 120°

To change the wedges, remove the padding from the brace. Place the appropriate wedge for the front and/or back in the side opening and screw the wedge in tight. Move the brace joint a few times to ensure that it can move freely and the wedges are secure. Mount the brace joint on the brace.

Only change the extension and flexion restrictions under the orders of the doctor treating you. The wedges set the exact radius of movement for the brace. If soft tissue conditions pose issues, then it may be necessary to use a larger wedge, due to soft tissue compression, in order to limit movement of the knee to the desired radius.

## Fitting instructions

- Loosen the tension of the cord load relief mechanism by pulling on the rotary seal in the upper section of the brace (fig. 1).
- Fold out the upper section of the brace using the silicone top band (fig. 2).
- Grip the brace with both hands on the sides at the upper end or grip the lateral donning devices (loops) on the inside (fig. 3). Bend your leg at approx. 70° and pull the brace onto the leg (fig. 4).
- Pull up the brace until the kneecap is positioned centrally within the thin comfort area. Please make sure that the brace has been placed on the leg straight and is not twisted.

- Now fold back the silicone top band in the upper section of the brace towards the skin side (fig. 6).
- Now tighten the cord load relief mechanism. First press in the rotary seal in the upper section (fig. 7). You will hear a clear click. A load-relieving cord tension builds up as a result of turning the wheel clockwise (fig. 8). Fully bend and stretch your knee two or three times and readjust the tension on the rotary seal if you need to. The cord tension loosens with millimetre precision as a result of turning the wheel anti-clockwise.

Depending on your level of pain, you can adjust the mechanism to fit your needs. Take a few steps and ensure that the brace is sitting correctly on the leg. If necessary, correct the position of the brace again.

### Taking the brace off

Loosen the tension of the cord load relief mechanism by pulling on the rotary seal in the upper section of the brace (Fig. 9). Fold out the silicone top band. Grip the brace with both hands on the lower end of the brace and pull it off the leg (fig. 10).

### Care instructions

Fats, oils, lotions, ointments and soap residues can affect the material and cause skin irritations.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

### Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

### Material composition

Polyester, polyamide, polyurethane, spring steel

### Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

### Disposal

Wishes you a speedy recovery!

### Your medi team

Wishes you a speedy recovery!

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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