**medi**

## lipomed bra

Kompressionsbekleidung  
compression garments

Gebrauchsanweisung. Instructions for use. Инструкция по использованию. Istruzioni per l'uso. Upute za uporabu. Bruksanvisning. دليل الاستخدام. Instrucciones de uso. Mode d'emploi. Gebruiksaanwijzing. הוראות שימוש. Naudojimo instrukcija. Instrukcja zakładania. 使用说明. Navodila za uporabo. Manual de instrucções. Käyttöohje. Інструкція з використання. Brugsvejledning. Instrucțiuni de utilizare.

**medi. I feel better.**

## Inhalt / Content

Deutsch / German	3–4
Englisch / English	4–5
Русский / Russian	5–6
Italiano / Italian	6–7
Hrvatski / Croatian	7–8
Norsk / Norwegian	8–9
عربي / Arabic	10–11
Español / Spanish	11–12
Français / French	12–13
Nederlands / Dutch	13–15
עברית / Hebrew	15–16
Lietuvių / Lithuanian	16–17
Polski / Polish	17–18
中文 / Chinese	18–19
Slovenščina / Slovenian	19–20
Português / Portuguese	20–21
Svenska / Swedish	22–23
Suomi / Finnish	23–24
Українська / Ukrainian	24–25
Dansk / Danish	25–26
Română / Romanian	26–27
Eesti / Estonia	27–28

# lipomed bra

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## Instructions for Use

### Intended purpose

Post-surgical garment for torso compression of tissue.

### Product features

The lipomed bra exerts even pressure on the tissue over the entire area. This compressive bra supports the breast and simultaneously reduces build-up of oedema and haematoma.

### Indications

All indications where compression in the torso area is necessary, such as:

- After breast augmentation
- After breast reduction
- After breast lift surgery
- After breast reconstruction

### Contraindications

- Arterial circulatory disorders
- Breathing difficulties
- Intolerance to material

### Risks and side effects

If swelling or paraesthesia are observed while wearing the compressive clothing, please contact your doctor immediately. You should only wear this product for long periods of time if recommended by your physician/Health Care Professional.

### Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that the persons supporting carer has received a briefing on this from healthcare professionals. Healthcare professionals should provide care to adults and children, applying the available information on measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

### Instructions for fitting and wear

The doctor or trained specialist must select the correct lipomed bra size based on the patient's measurements. If the patient has undergone breast augmentation, it is recommended that this product be used in conjunction with the lipomed belt in order to prevent implant displacement. Lipomed products are designed for unique use by one patient (lasting 4 - 6 weeks) following an operation (single patient use). If the product is used in multiple instances or if it is modified, then its medical safety is no longer guaranteed. Changes to the design, cut and standard model can result in impairments in the therapeutic medical efficacy of the product and may put the patient at risk. Changes of any kind executed by third parties shall render the EU declaration of conformity (CE marking) required by law for medical devices which is issued by the manufacturer null and void.

Contact with chemical substances (e.g. benzenes, methanols, salicylic acids in high concentrations) may impair the compressive function of the products. Oils and ointments should not be applied to when compression garments are used.

### Length of time that compression clothing should be worn and used for

The doctor will stipulate the period of time that the device needs to be worn for depending on the type and extent of the surgery. Generally, the device must be worn for around four weeks. To elicit optimal results, the lipomed bra should be worn day and night.

### Instructions for care

The lipomed compression bra can be washed in a normal washing machine, using a mild, non-bleaching detergent. The wash temperature should not exceed 40°C. Do not use fabric conditioner. Wash your lipomed compression clothing separately or with clothing that is similar in colour.

**Note:** Before every wash, close any hook-and-loop fasteners in order to preserve their functionality and to prevent damage to other items of clothing. Air dry, but

not in direct sunlight or on heating elements. Excessively high temperatures may damage the elasticity of the material and thus the medical efficacy of the lipomed bra. In order to prevent longer interruptions in compression therapy, it is recommended that the patient procures another lipomed bra in the correct size so that they can change between the two bras.

### Storage

The products must be stored in a dry place and out of direct sunlight.

### Composition:

51% polyamide  
35% cotton  
14% elastane

### Disposal

Please dispose of your compression wear via your general waste or household waste. There are no specific criteria for disposing of the compression wear.

### Complaints

In the event of any complaints regarding the product such as damage to the knitted fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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## Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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medi UK Ltd.