medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-57
info@medi.de
www.medi.de

medi Australia Pty Ltd
 Unit 4/13-21 Hallmark Street
 Pendle Hill NSW 2145
 Australia
 T +61 2 9890 8696
 F +61 2 9890 8439
 sales@medi.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 verttrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: +32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Canada Inc / médi Canada Inc
 104-1375 Lionel-Boulet,
 Varennes, Québec,
 QC Canada J3X 1P7
 T +1 450-583-3317 / +1 800 361 3153
 F +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

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embraceME / Bratelle

Intended purpose

Post-surgical garment for torso compression  of tissue.

Indications

All indications where compression in the torso area is necessary, such as:

- Edema in the torso area / upper body
- Edema following breast surgery
- Scar tissue

Contraindications

- advanced peripheral arterial occlusive disease
- Decompensated heart failure (NYHA III + IV)
- Septic phlebitis
- Phlegmasia coerulea dolens

Cautions

- Scar tissue after laparotomy or thoracotomy
- Untreated, local or systemic infections
- Respiratory problems
- Open wounds in the compression area
- Pronounced weeping dermatoses
- Intolerance to material(s)
- Severe sensitivity disorders
- Advanced peripheral neuropathy (e.g. as seen diabetes mellitus)
- Primary chronic polyarthritis

Important note

The correct measurement for embraceME / Bratelle must be carried out by a doctor or specialist.

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for

this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Instructions for putting embraceME on

Fig. 1 shows your embraceME from the inside. There you will find a sewn-in pocket. The arrow in the illustration points to the opening of the pocket into which the molding is pushed.

Fig. 2 shows the molding. There is an arrow and two letters on it. The arrow shows the insertion direction of the molding. The letters indicate the size of the molding and the right (R) or left (L) side. In this example, the M in "M - R" stands for size M and the R for the right side.

Fig. 3: Push the molding into the opening with the arrow marking first. As soon as it is completely in the pocket, unfold the molding again and, if its necessary, reposition it so that it sits comfortably in the pocket.

Fig. 4 shows you embraceME from the inside with the inserted moldings. The markings make it clear again how to insert it correctly.

Fig. 5 Put on your embraceME. Please ensure that the horizontal seam of the product sits directly below the bust for a proper fit. Close the hook-and-loop closure starting at the hip and moving upward toward the neck.

Fig.6 Hold the side of the zipper with the slider firmly in one hand to keep the slider stable and prevent it from twisting. Insert the flat, rigid end of the other zipper half straight and fully into the opening at the bottom of the slider. Ensure that the pin is pushed in all the way and that both fabric edges align neatly. Keep holding the lower section firmly and pull the slider upward in a smooth motion until the zipper is completely closed.

Fig.7 Position the zipper cover so that it lies flat, does not shift, and fully covers the zipper.

Fig.8 Wear your clothing over embraceME. The product achieves its full effectiveness only with direct skin contact. Therefore, clothing must not be worn between the skin and the product.

Instructions for putting Bratelle on

Fig.1 Bring the two sides of the product together at the front and close the hook-and-loop closure. Make sure the closure noticeably clicks into place.

Fig.2 Rotate the product toward the spine and move it to the back. The closure should be positioned in the center of the spine.

Fig.3 Align the cups so that they fully enclose each breast. Make sure the border band lies horizontally and smoothly against the body. If needed, pull the border band downward so it sits securely beneath the bust.

Fig.4 Pull both straps over your shoulders.

Fig.5 Check whether the straps are neither pinch nor slipping, and adjust them as needed.

Fig.6 The product should now fit comfortably, provide support, and neither press nor create wrinkles.

Instructions for care

Wash your embraceME / Bratelle Compression breast bandage compression clothing separately or with clothing that is similar in colour. We recommend using a laundry bag. Before washing, remove the moldings from the embraceME to prevent damage. Fasten the hook-and-eye closure before each wash to protect the material and other clothing from damage. Air-dry, but not in direct sunlight or on heating elements. Excessively high temperatures may damage the elasticity material and thus the medical efficacy of the embraceME / Bratelle. In order to prevent longer interruptions in compression therapy, it is recommended that the patient procures another embraceME / Bratelle so that they can change between the two shirts.

-  Wash the product by hand, preferably with medi clean soap or in a delicates wash cycle at 40°C with mild detergent and without fabric conditioner.
-  Do not bleach.
-  Air dry or dry in a tumble dryer using the delicates cycle.
-  Do not iron.
-  Do not dry-clean.

Storage

The products must be stored in a dry place and out of direct sunlight.

Length of time that the products should be worn and used for

6 months – due to material wear and tear, medical effectiveness can only be

guaranteed for a defined period of use.
The requires correct handling (e.g. when caring for, putting on and taking off).

Disposal

Dispose of properly after use.

Liability

The manufacturer's liability will become void if the product is not used as intended (including own repairs or those of external service providers).

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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Important notes

This medical device **[MD]** is made only for multiple uses on a single patient . If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after thorough consultation with a medical professional.

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