medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
Unit 4/13-21 Hallmark Street
Pendle Hill NSW 2145
Australia
T +61 2 9890 8696
F +61 2 9890 8439
sales@medi.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43 512 57 95 15
F +43 512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: + 32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Canada Inc / médi Canada Inc
104-1375 Lionel-Boulet,
Varennes, Québec,
QC Canada J3X 1P7
T +1 450-583-3317 / +1 800 361 3153
F +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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
medi Trading (Shanghai) Co. Ltd.
Room 613, HnR Group Building
No.1178, Tian Yao Qiao Road
Xuhui District, Shanghai, 200232
People's Republic of China
T +86-21 50582318
infocn@medi.cn
www.medi.cn

UK Responsible Person
medi UK Ltd.
Plough Lane
Hereford HR4 OEL
Great Britain
T +44-1432 37 35 00
F +44-1432 37 35 10
enquiries@mediuk.co.uk
www.mediuk.co.uk

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Important notes
This medical device MD is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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English

medi 4C flex

Intended purpose
The medi 4C flex is a brace for relieving load on the lumbar spine/thoracic spine (sagittal and frontal limitation of range of motion)

Indications
All indications requiring a partial relief by correction of the lumbar/thoracic spine and movement restriction in sagittal and frontal plane, such as:
- Stable vertebral compression fractures of the middle and lower thoracic and lumbar spine
- Temporary postoperative immobilisation
- Nonsurgical inoperable or postoperative treatment for tumours and metastases of the spinal column
- Axial stabilisation for osteoporotic vertebral collapse in the thoracic and lumbar spine

Contraindications
Unstable fractures of the thoracic / lumbar spine with and without neurological deficit

Risks / Side Effects
If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:
- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

Intended users and patient target groups
Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.
Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/ indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Fitting instructions
Please wear the product on thin outerwear.
- **Adjusting the width**
 First loosen the screws on the chassis and adjust the orthosis to the width of the patient's body. Please always check that the functioning of the hinged fastener is not restricted.
 Now loosen the two screws on the side of the chassis and the screws on the fastener. Adjust the orthosis to suit the patient's length – distance from the symphysis pubis to the subclavicular contact area. Then retighten all the screws.
- **Adjusting the strap length**
 Using the hook and loop fastener and the reversing loop adjust the length of the strap so that the pad comes to lie in the middle of the back. Remember that a few centimetres of the strap will have to be pulled through the fastener.
- **Opening and closing the fastener**
 Slide the grey bar sideways to open the fastener. To close the orthosis flip the

cover over until it audibly clicks into place.

Care instructions
The orthosis can be cleaned with a damp cloth, the back pad with a brush. Soap residues, lotions and ointments can cause skin irritation and material wear.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Storage instructions
Keep the product in a dry place and do not expose to direct sunlight.

Material composition
Aluminium, PVC

Liability
The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal
Dispose of properly after use.

Your medi Team
Wishes you a quick return to full fitness.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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