medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
www.medi.de

medi Australia Pty Ltd
Unit 4/13-21 Hallmark Street
Pendle Hill NSW 2145
Australia
T +61 2 9890 8696
F +61 2 9890 8439
sales@medi.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43 512 57 95 15
F +43 512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: + 32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Canada Inc / médi Canada Inc
104-1375 Lionel-Boulet,
Varennes, Québec,
QC Canada J3X 1P7
T +1 450-583-3317 / +1 800 361 3153
F +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.
Room 613, HnR Group Building
No.1178, Tian Yao Qiao Road
Xuhui District, Shanghai, 200232
People's Republic of China
T +86-21 50582318
infocn@medi.cn
www.medi.cn

UK Responsible Person
medi UK Ltd.
Plough Lane
Hereford HR4 OEL
Great Britain
T +44-1432 37 35 00
F +44-1432 37 35 10
enquiries@mediuk.co.uk
www.mediuk.co.uk

4 049772 425945

E011935/07.2025

circaid® reduction kit™ lower leg inelastic compression system

Gebrauchsanweisung. Instructions for use. Mode d'emploi. Instrucciones de uso. Instruções para aplicação. Istruzioni per l'uso. Gebruiksaanwijzing. Bruksanvisning. Инструкция по использованию. Instrukcja zakładania. دليل الاستخدام. Használati útmutató. 使用说明. הוראות שימוש. Bruksanvisning. Kasutusjuhend. Käyttöohje. 使用説明書.

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Important notes
This medical device MD is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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circaid reduction kit lower leg

Intended purpose
The compression system is designed to provide compression to the leg for patients with venous and lymphatic disorders.

Indications
- Lymphedema
- Lipedema
- Chronic venous insufficiency
 - Varicose veins
 - Lipodermatosclerosis
 - Venous leg ulcer
- Deep vein thrombosis/Thrombosis Prevention
- Acute deep vein thrombosis
- Post Thrombotic Syndrome
- Post Sclerotherapy
- Dependent oedema

Contraindications
- Advanced peripheral arterial disease (if one of these parameters is present: ABPI < 0.5; ankle arterial pressure < 60mmHg; toe pressure < 30 mmHg or TcPO2 < 20 mmHg on dorsum of foot). If non-elastic materials are used, the patient can still try out compression clothing if they have ankle arterial pressure between 50 and 60 mmHg, under close clinical supervision
- Decompensated heart failure
- Septic phlebitis
- Phlegmasia cerulea dolens
- Untreated infection of the leg and/or foot
- Any circumstances under which an increase in venous or lymphatic reflux is undesirable
- Suspected or known untreated acute deep vein thrombosis

Risks / Side Effects
- Mild to Moderate Peripheral Arterial Disease
- Impaired Sensation - must be able to sense applied pressure
- Intolerance to material(s)
- Pediatric & assisted use - must be able to sense & communicate applied pressure

Intended users and patient target groups
Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.
Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Individualisation
- Separate the body and spine from one another, unroll the bottom band and locate the number markings along the bottom and top edge of the body piece.
- Locate the patient's ankle circumference measurement amongst the printed numbers along the bottom edge. Then, locate the patient's calf circumference measurement amongst the printed numbers along the top edge of the garment.
- Take the spine and place the bottom left edge on the appropriate number. Take the top top left edge of the spine and place it on the appropriate number. Adjust the center spine tab to create a smooth spine curve.

Note: *There is no one correct spine shape. The garment may not lay flat once the spine has been set. This is to be expected when dealing with more extreme spine angles. If there is a large difference between the two circumferences, the spine may not reach both the bottom and top edges of the body piece. In this case, center the spine. When dealing with more extreme spine angles, excess material sticking out beyond the top of the spine piece can be trimmed.*

- With the spine in place, fold over the spine bands exposing the body piece. Cut the excess material from the body piece leaving approximately 5 cm of overlap. Be sure not to cut any other part of the garment.

Note: *If desired, attach the separate hook tab stays along the spine for additional adhesion. You may also trace the spine to easily reattach it in its original position in the unlikely event of the spine becoming detached.*

- If the lower leg component is too long for the patient, it can be trimmed by 3 cm. Trim along the indicated cut line on the black side of the material. Trim both the body and the spine portion.

Donning instructions
- If the foot or toes require any bandaging, apply prior to donning the undersleeve and legging.
- Slide the enclosed circaid undersleeve or undersock onto the leg up to the knee (for lower leggings) or to the groin (for upper leg and whole leg). The undersleeve may be trimmed in length if desired. If a knee and/or upper leg component will be applied it is not recommended to trim the undersleeve until these components have been fitted over the undersleeve to determine the excess length.
- Align the garment just above the ankle bone and below the knee crease. The black side of the material should face the skin. Unroll the top two bands and loosely secure them to the outside of the garment to anchor it in place. This enables you to easily raise or lower the garment for correct positioning.
- Unroll the bottom two bands and secure the bottom most band to the outside of the garment to a firm and comfortable compression level.
- While holding the second band, detach the next band. Secure the second band and continue this process with the other bands, going up the leg.
- Ensure that the bands are placed in an alternating order.
- Once all bands are secured, inspect the garment for any gaps or creases. Adjust the bands as necessary.
 If parts of the leg require additional support or a gap exists, use the provided shelf strap as needed.

Note: *If no reduction kit knee component or upper leg component are being worn, you may fold the top and bottom part of the circaid undersleeve over the top and bottom edge of the circaid compression garment.*

- Apply an appropriate foot compression garment as needed. If using another circaid garment, please refer to the specific product instructions included in its packaging.

Using the Built-In-Tension System (BTS)
Bands should be pulled to a firm but comfortable tension to provide therapeutic compression within the tolerance level of the patient.
- Starting with the bottom band, line up the black triangle on the card with one of the BTS lines on the bottom band.
- Note where the second BTS line lines up with the card´s compression ranges.
- If the BTS line on the garment either falls short or goes beyond the desired compression range, readjust the band as necessary so that the second BTS line on the garment is aligned with the desired compression range. Completely loosen the band before reapplying it if too much tension has been applied.

Note: *Record tension setting for each set of markings so the patient can easily repeat the process and ensure it is donned to the same or higher tension.*

- Repeat the steps for each band going up the leg.
 It is not required to have all bands set to the same tension level. If more than one level of tension is appropriate for treatment, use the BTS card to identify and record the appropriate tensions levels.

Bands may need adjusting throughout the day to maintain the desired compression range as edema is reduced.

Care instructions
Hook and loop fasteners should be closed for washing. This compressive product can be washed on a delicate cycle and dried in a tumble dryer set to a low heat. We recommend using a mesh laundry bag. To extend the life of the product, hand wash and drip dry.

Machine wash warm - gentle cycle
Do not bleach
Tumble dry, low heat
Do not iron
Do not dry clean
Do not wring

Storage instructions
Keep the product in a dry place and do not expose to direct sunlight.

Lifespan
2 months – due to material wear and/or slackening, medical efficacy can only be guaranteed for a defined lifespan. This assumes correct handling of the product (e.g. during care and when putting on and taking off the product).

Material composition
Polyurethane, Nylon, Elastane

Disposal
The product can be disposed of in the domestic waste.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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