medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth / Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
Unit 4/13-21 Hallmark Street
Pendle Hill NSW 2145
Australia
T +61 2 9890 8696
F +61 2 9890 8439
sales@medi.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43 512 57 95 15
F +43 512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: + 32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Canada Inc / médi Canada Inc
104-1375 Lionel-Boulet,
Varennes, Québec,
QC Canada J3X 1P7
T +1 450-583-3317 / +1 800 361 3153
F +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.
Room 613, HnR Group Building
No.1178, Tian Yao Qiao Road
Xuhui District, Shanghai, 200232
People's Republic of China
T +86-21 50582318
infocn@medi.cn
www.medi.cn

Medi Turk Ortopedi Medikal Ic ve Dis Tic. AS
Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5
06510 Çankaya Ankara
Turkey
T: +90 312 435 20 26
F: +90 312 434 17 67
info@medi-turk.com
www.medi-turk.com

medi Ukraine LLC
Evhena Sverstiuka str, 11.
Kiev 02002
Ukraine
T: +380 44 591 11 63
F: +380 44 392 73 73
info@medi.ua
http://medi.ua

medi USA L.P.
6481 Franz Warner Parkway
Whitsett, N.C. 27377-3000
USA
T +1-336 4 49 44 40
F +1-888 5 70 45 54
info@mediusa.com
www.mediusa.com

UK Responsible Person
medi UK Ltd.
Plough Lane
Hereford HR4 OEL
Great Britain
T +44-1432 37 35 00
F +44-1432 37 35 10
enquiries@mediuk.co.uk
www.mediuk.co.uk

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Important notes
This medical device MD is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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English

Stabimed® RICE

Intended purpose
The Stabimed RICE is a knee support brace with limitation of extension / flexion.

Indications
All indications that require reinforcement of the physiological stabilisation of the knee joint and/or stress relief of the joint, such as:
- Following injury to the collateral ligament
- Following meniscal injury
- Following endoprostehtic intervention
- For knee instability
- For knee osteoarthritis

Contraindications
None known at present.

Risks / Side Effects
If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:
- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material. The immobilisation which is achievable by the device presents a risk for development of thrombosis in the lower veins of the legs, especially in the context of surgical interventions or trauma.

Intended users and patient target groups
Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.
Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/ indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Setting the limits for flexion and extension (can only be carried out by an orthopaedic technician)
Extension restriction wedges, E (silver): (stretch restriction): 0°, 10°, 20°, 30°, 45°
Flexion restriction wedges, F (black) (bend restriction): 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 120°
The medial (internal) and lateral (external) degree settings must always be the same.
1. If you do not install the wedges in the knee brace, it will automatically be set to 0° extension and flexion is not limited.
2. To change the wedges, you need to turn round the joint so that the joint cover is pointing downwards. Select the desired limitation using the wedges included in the delivery. The relevant degree setting is stamped on the wedges.
3. Once you have inserted the wedge, secure it by screwing in the screw using the included screwdriver.
4. Move the brace joints and check that

the wedges are securely in place.

Important notes
Only change the extension and flexion limits on instruction from your doctor or orthotist.

Fitting instructions
- Open both straps (Fig. 1). The valve is at the top right and the straps are also facing right.
- Bend your leg to approximately 45°. Place the brace onto the front of the leg so that the lower edge of the hole for the kneecap is sitting on the lower side of the kneecap and is positioned centrally (Fig. 2).
- Now position both of the joints on the brace. To do this, put the joint splints in the small attached hook-and-loop pockets (Fig. 3).
- Now stick the splints onto the side of the brace. Ensure that the brace joints are in the correct position: The middle of the brace joints should be level with the top of the kneecap (Fig. 4). The joints must be positioned behind the centre of the side of the leg and must run parallel to each other. The joints must be adapted to the patient's individual leg shape by limiting them.
- Wrap both straps around the leg and close the straps on the brace (Fig. 5).
- If the straps are too short, simply fasten both of the additional strap extensions on the hook-and-loop areas at the ends of the straps for the brace (Fig. 6). You can now wrap the straps around the leg as usual and fasten them on the brace.

Care instructions
Please close the hook and loop fasteners before washing and remove the splint. Switch the valve to **"OFF"**. We recommend using a laundry bag. Soap residues, lotions and ointments can cause skin irritation and material wear. The cooling compress can be cleaned carefully with water, a moistened towel or a sponge.
- Wash the product by hand, preferably using a medi clean detergent, or in delicate cycle at 30°C using a mild detergent without fabric conditioners.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Storage instructions
Keep the product in a dry place and do not expose to direct sunlight.

Material composition
Knee sleeve: Cotton, nylon
Blow-up inner cushion: Polyester, nylon
Air pump: Silicone, polyethene, synthetic rubber
Cooling compress: Propylene glycol, sodium carboxymethyl cellulose, water
Joints: Aluminium, polyester

Liability
The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal
Dispose of properly after use.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section

2, No. 65 of the Regulation (EU) 2017/745 (MDR).

Note on additional functions of the brace
The brace can also be used with a cooling compress and/or with compression with the aid of an additional air pump, dependent on the instructions of the doctor.

Warning!
Please consult the doctor treating you before applying compression and/or cold therapy. This is mandatory. It may be that this treatment is contraindicated for you.
If the cooling compress is defective for any reason, if it leaks and the gel comes into contact with the skin or eyes, wash the affected areas thoroughly with water straight away. Never try to heat the cooling compress in any way (e.g. water baths, microwaving, ovens etc. are all prohibited).

Contraindications for cooling and compression therapy
If you have a vascular disease, circulation disorder, paraesthesia or tissue inflammation, these therapies may only be introduced once they have been discussed with your doctor.

Additional function: Cooling
- To use the cooling function, leave the included cooling compress in the freezer briefly (max. 30 minutes) or simply in the fridge. An excessively long cooling period in the fridge is not recommended as the temperature could then drop too low and injuries or frostbite could occur.
- After successful cooling, the cooling compress (Fig. 7) can simply be fastened to the black inner side of the brace (the blue inner side of the cooling compress is left facing you).

Warning!
Ensure that an underlay (thin towel or piece of clothing) is always present between the skin and the cooling compress. The length of time and frequency of the cooling phases must always be specified by the treating doctor and should not last longer than a max. of 30 minutes.

Additional function: Compression
- To apply compression, attach the tube to the valve (upper right) (Fig. 8).
- Switch the valve to "**ON**" (Fig. 9).
- Pump carefully until you have achieved the desired level of compression (Fig. 10). After filling has been completed, switch the valve to "**OFF**".
- To let some air out, you can press the button on the pump (Fig. 11).
- If you want to let all of the air out, remove the tube and set the valve setting to "**ON**". Now the air can escape again.

Warning!
Do not overfill the air chamber.
Excessive pressure can cut off the blood vessels and can lead to decreased circulation of blood and other fluids.

Your medi Team
Wishes you a speedy recovery!

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