medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth / Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de/en

medi Australia Pty Ltd
Unit 4/13-21 Hallmark Street
Pendle Hill NSW 2145
Australia
T +61 2 9890 8696
F +61 2 9890 8439
sales@medi.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43 512 57 95 15
F +43 512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: + 32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Canada Inc / médi Canada Inc
104-1375 Lionel-Boulet,
Varennes, Québec,
QC Canada J3X 1P7
T +1 450-583-3317 / +1 800 361 3153
F +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.
Room 613, HnR Group Building
No.1178, Tian Yao Qiao Road
Xuhui District, Shanghai, 200232
People's Republic of China
T +86-21 50582318
infocn@medi.cn
www.medi.cn

UK Responsible Person
medi UK Ltd.
Plough Lane
Hereford HR4 OEL
Great Britain
T +44-1432 37 35 00
F +44-1432 37 35 10
enquiries@mediuk.co.uk
www.mediuk.co.uk

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Important notes
This medical device MD is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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Polyamid, Polyester, Elastan, Viskose

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Levamed active

Intended purpose
Levamed active is a support for compression of the soft tissue in the ankle joint with additional functional elements.

Indications
All indications in which soft tissue compression with local/local additional pad(s) at the ankle joint is necessary and at least one direction of movement is restricted by additional functional elements, such as:
- Joint sprain with overstretch of the capsule
- Chronic ligamentous instability with recurrent supination traumas (twisting of the foot)
- In the second phase of treatment for lateral ligament injuries
- Nonsurgical treatment of injuries of the syndesmosis not requiring surgery

Contraindications
None known at present.

Risks / Side Effects
If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:
- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

Intended users and patient target groups
Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.
Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Note on wearing
The support achieves its greatest effect during physical activity. In general, the brace can be worn all day. However, this should be done according to your comfort. During longer breaks (e.g. prolonged sitting / road trips, sleeping), the brace should be removed.

Fitting instructions
- Take hold of the top of the support on both sides and pull it over your foot until the heel part is sitting correctly and the silicone inserts are positioned

around the ankle. (Fig. 1)
- Place the attached strap on the floor in front of you, with the side with the footprint on it facing up. (Fig. 2)
- Now position your foot on the symbol so that the strap runs under the middle of the sole of your foot. Ensure that your foot is at 90° to your shin, as far as is possible. (Fig. 3)
- Take hold of the end of the strap, which is on the inside of your foot ◰, and wrap it over the upper side of your foot diagonally so that it is now on the outer side of your foot. Fasten the end of the strap above the outer side of your ankle, using the hook and loop area intended for this purpose. (Fig. 4)
- Then take the shorter strap on the outside of your foot ◰ and fasten it above the outer side of your ankle too. (Fig. 5) Pull on this end of the strap to secure the entire belt around your foot. This ensures comfortable stability for the joint.
- Finally, take the longer strap ◰ on the outside of your foot across the top of your foot diagonally and then around the leg, bringing it to the front. Stick the end of the strap to the marking on the strap under it. (Fig. 6+7+8)

NOTE: If you need to, you can shorten the strap to fit you by cutting it at the markings ◰ or ◰. To do this, unfasten the hook and loop fasteners at the end of the strap and shorten it accordingly.

Care instructions

Hook and loop fasteners should be closed for washing. Fabric softeners, fats, oils, lotions, ointments and soap residues can affect the material, hinder the Clima Comfort effect and cause skin irritations.
- Wash the product by hand, preferably using a medi clean detergent, or in delicate cycle at 30°C using a mild detergent without fabric conditioners.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Material composition

Polyamide, polyester, elastane, viscose

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

Dispose of properly after use.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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